

CM LOC® Kompatibilitätsliste / list of compatible system and available abutments

Implant(at) System	Durchmesser / Diameters	Wenn benötigt, Implantat-Adaptor As needed, implant adapter	Bestell-Nr. / order n°	Torque abutment
Astra tech OsseoSpeed®	D 3.5, 4.0 mm	CM LOC® abutment OsseoSpeed® 3.5/4.0 H1 CM LOC® abutment OsseoSpeed® 3.5/4.0 H2 CM LOC® abutment OsseoSpeed® 3.5/4.0 H3 CM LOC® abutment OsseoSpeed® 3.5/4.0 H4 CM LOC® abutment OsseoSpeed® 3.5/4.0 H5 CM LOC® case guide OsseoSpeed® 3.5/4.0 20°	05001452 05001453 05001454 05001455 05001456 05001587	25 Ncm -
	D 4.5, 5.0 mm	CM LOC® abutment OsseoSpeed® 4.5/5.0 H1 CM LOC® abutment OsseoSpeed® 4.5/5.0 H2 CM LOC® abutment OsseoSpeed® 4.5/5.0 H3 CM LOC® abutment OsseoSpeed® 4.5/5.0 H4 CM LOC® abutment OsseoSpeed® 4.5/5.0 H5 CM LOC® case guide OsseoSpeed® 4.5/5.0 20° CM LOC® FLEX abutment OsseoSpeed® 4.5/5.0 H1 CM LOC® FLEX abutment OsseoSpeed® 4.5/5.0 H2 CM LOC® FLEX abutment OsseoSpeed® 4.5/5.0 H3 CM LOC® FLEX abutment OsseoSpeed® 4.5/5.0 H4 CM LOC® FLEX abutment OsseoSpeed® 4.5/5.0 H5 CM LOC® FLEX case guide OsseoSpeed® 4.5/5.0 30°	05001457 05001458 05001459 05001460 05001461 05001588 05001770 05001771 05001772 05001773 05001774 05001894	- -
Astra tech OsseoSpeed® EV	D 3.6 mm	CM LOC® abutment OsseoSpeed EV 3.6 H1 CM LOC® abutment OsseoSpeed EV 3.6 H2 CM LOC® abutment OsseoSpeed EV 3.6 H3 CM LOC® abutment OsseoSpeed EV 3.6 H4 CM LOC® abutment OsseoSpeed EV 3.6 H5 CM LOC® case guide OsseoSpeed EV 3.6 20°	05001963 05001964 05001965 05001966 05001967 05001981	25 Ncm -
	D 4.2 mm	CM LOC® abutment OsseoSpeed® EV 4.2 H1 CM LOC® abutment OsseoSpeed® EV 4.2 H2 CM LOC® abutment OsseoSpeed® EV 4.2 H3 CM LOC® abutment OsseoSpeed® EV 4.2 H4 CM LOC® abutment OsseoSpeed® EV 4.2 H5 CM LOC® case guide OsseoSpeed® EV 4.2 20° CM LOC® FLEX abutment OsseoSpeed® EV 4.2 H1 CM LOC® FLEX abutment OsseoSpeed® EV 4.2 H2 CM LOC® FLEX abutment OsseoSpeed® EV 4.2 H3 CM LOC® FLEX abutment OsseoSpeed® EV 4.2 H4 CM LOC® FLEX abutment OsseoSpeed® EV 4.2 H5 CM LOC® FLEX case guide OsseoSpeed® EV 4.2 20°	05991968 05001969 05001970 05001971 05001972 05001982 05001932 05001933 05001934 05001935 05001936 05001946	- -
	D 4.8 mm	CM LOC® abutment OsseoSpeed® EV 4.8 H1 CM LOC® abutment OsseoSpeed® EV 4.8 H2 CM LOC® abutment OsseoSpeed® EV 4.8 H3 CM LOC® abutment OsseoSpeed® EV 4.8 H4 CM LOC® abutment OsseoSpeed® EV 4.8 H5 CM LOC® case guide OsseoSpeed® EV 4.8 20° CM LOC® FLEX abutment OsseoSpeed® EV 4.8 H1 CM LOC® FLEX abutment OsseoSpeed® EV 4.8 H2 CM LOC® FLEX abutment OsseoSpeed® EV 4.8 H3 CM LOC® FLEX abutment OsseoSpeed® EV 4.8 H4 CM LOC® FLEX abutment OsseoSpeed® EV 4.8 H5 CM LOC® FLEX case guide OsseoSpeed® EV 4.8 20°	05001973 05001974 05001975 05001976 05001977 05001983 05001937 05001938 05001939 05001940 05001941 05001947	- -
AUROSAN REVOIS® PRO	D 3.8, 4.3, 5.0 mm	CM LOC® abutment REVOIS® PRO H1 CM LOC® abutment REVOIS® PRO H2 CM LOC® abutment REVOIS® PRO H3 CM LOC® abutment REVOIS® PRO H4 CM LOC® abutment REVOIS® PRO H5 CM LOC® case guide REVOIS® PRO 20° CM LOC® FLEX abutment REVOIS® PRO H1 CM LOC® FLEX abutment REVOIS® PRO H2 CM LOC® FLEX abutment REVOIS® PRO H3 CM LOC® FLEX abutment REVOIS® PRO H4 CM LOC® FLEX abutment REVOIS® PRO H5 CM LOC® FLEX case guide REVOIS® PRO 30°	MD6001871 MD6001872 MD6001873 MD6001874 MD6001875 MD6001876 MD6001881 MD6001882 MD6001883 MD6001884 MD6001885 MD6001886	25 Ncm -

Camlog® Camlog®	D 3.8 mm	CM LOC® abutment Camlog® 3.8 H1 CM LOC® abutment Camlog® 3.8 H2 CM LOC® abutment Camlog® 3.8 H3 CM LOC® abutment Camlog® 3.8 H4 CM LOC® case guide Camlog® 3.8 20°	05001544 05001545 05001546 05001547 05001591	30 Ncm -
	D 4.3 mm	CM LOC® FLEX abutment Camlog® 3.8 H1 CM LOC® FLEX abutment Camlog® 3.8 H2 CM LOC® FLEX abutment Camlog® 3.8 H3 CM LOC® FLEX abutment Camlog® 3.8 H4 CM LOC® FLEX case guide Camlog® 3.8 30° CM LOC® abutment Camlog® 4.3 H1 CM LOC® abutment Camlog® 4.3 H2 CM LOC® abutment Camlog® 4.3 H3 CM LOC® abutment Camlog® 4.3 H4 CM LOC® case guide Camlog® 4.3 20° CM LOC® FLEX abutment Camlog® 4.3 H1 CM LOC® FLEX abutment Camlog® 4.3 H2 CM LOC® FLEX abutment Camlog® 4.3 H3 CM LOC® FLEX abutment Camlog® 4.3 H4 CM LOC® FLEX case guide Camlog® 4.3 30°	05001790 05001791 05001792 05001793 05001897 05001549 05001550 05001551 05001552 05001592 05001795 05001796 05001797 05001798 05001898	- -
Camlog® Conelog®	Conelog® D 3.8, 4.3 mm	CM LOC® abutment Conelog® H1 CM LOC® abutment Conelog® H2 CM LOC® abutment Conelog® H3 CM LOC® abutment Conelog® H4 CM LOC® abutment Conelog® H5 CM LOC® case guide Conelog® 20° CM LOC® FLEX abutment Conelog® H1 CM LOC® FLEX abutment Conelog® H2 CM LOC® FLEX abutment Conelog® H3 CM LOC® FLEX abutment Conelog® H4 CM LOC® FLEX abutment Conelog® H5 CM LOC® FLEX case guide Conelog® 30°	05001482 05001483 05001484 05001485 05001486 05001594 05001805 05001806 05001807 05001808 05001809 05001899	30 Ncm -
Dentsply Ankylos® C	D 3.5, 4.5, 5.5, 7.0 mm	CM LOC® abutment Ankylos® C H1 CM LOC® abutment Ankylos® C H2 CM LOC® abutment Ankylos® C H3 CM LOC® abutment Ankylos® C H4 CM LOC® abutment Ankylos® C H5 CM LOC® case guide Ankylos® C 20° CM LOC® FLEX abutment Ankylos® C H1 CM LOC® FLEX abutment Ankylos® C H2 CM LOC® FLEX abutment Ankylos® C H3 CM LOC® FLEX abutment Ankylos® C H4 CM LOC® FLEX abutment Ankylos® C H5 CM LOC® FLEX case guide Ankylos® C 30°	05002014 05002015 05002016 05002017 05002018 05002020 05002002 05002003 05002004 05002005 05002006 05002008	25 Ncm -
Nobel Biocare Replace Select NP	D 3.5 mm	CM LOC® abutment Replace Select® NP3.5 H1 CM LOC® abutment Replace Select® NP3.5 H2 CM LOC® abutment Replace Select® NP3.5 H3 CM LOC® abutment Replace Select® NP3.5 H4 CM LOC® abutment Replace Select® NP3.5 H5 CM LOC® case guide Replace Select® NP3.5 20°	05001423 05001424 05001425 05001426 05001427 05001581	35 Ncm -
Nobel Biocare Replace Select RP	D 4.3 mm	CM LOC® abutment Replace Select® RP4.3 H1 CM LOC® abutment Replace Select® RP4.3 H2 CM LOC® abutment Replace Select® RP4.3 H3 CM LOC® abutment Replace Select® RP4.3 H4 CM LOC® abutment Replace Select® RP4.3 H5 CM LOC® case guide Replace Select® RP4.3 20° CM LOC® FLEX abutment Replace Select® RP4.3 H1 CM LOC® FLEX abutment Replace Select® RP4.3 H2 CM LOC® FLEX abutment Replace Select® RP4.3 H3 CM LOC® FLEX abutment Replace Select® RP4.3 H4 CM LOC® FLEX abutment Replace Select® RP4.3 H5 CM LOC® FLEX case guide Replace Select® RP4.3 30°	05001296 05001300 05001301 05001302 05001303 05001582 05001740 05001741 05001742 05001743 05001744 05001888	35 Ncm -
Nobel Biocare Active NP	D. 3.5 mm	CM LOC® abutment Active® NP3.5 H1 CM LOC® abutment Active® NP3.5 H2 CM LOC® abutment Active® NP3.5 H3 CM LOC® abutment Active® NP3.5 H4	05002109 05002110 05002111 05002112	35 Ncm

		CM LOC® abutment Active® NP3.5 H5 CM LOC® case guide Active® NP3.5 20°	05002113 05002187	-
Nobel Biocare Active RP	D 4.3 mm	CM LOC® abutment Active® RP4.3 H1 CM LOC® abutment Active® RP4.3 H2 CM LOC® abutment Active® RP4.3 H3 CM LOC® abutment Active® RP4.3 H4 CM LOC® abutment Active® RP4.3 H5 CM LOC® case guide Active® RP4.3 20° CM LOC® FLEX abutment Active® RP4.3 H1 CM LOC® FLEX abutment Active® RP4.3 H2 CM LOC® FLEX abutment Active® RP4.3 H3 CM LOC® FLEX abutment Active® RP4.3 H4 CM LOC® FLEX abutment Active® RP4.3 H5 CM LOC® FLEX case guide Active® RP4.3 30°	05001437 05001438 05001439 05001440 05001441 05001584 05001750 05001751 05001752 05001753 05001754 05001890	35 Ncm -
Nobel Biocare Brånemark® RP	D 4.0 mm	CM LOC® abutment Brånemark® RP4.0 H1 CM LOC® abutment Brånemark® RP4.0 H2 CM LOC® abutment Brånemark® RP4.0 H3 CM LOC® abutment Brånemark® RP4.0 H4 CM LOC® abutment Brånemark® RP4.0 H5 CM LOC® case guide Brånemark® RP4.0 20° CM LOC® FLEX abutment Brånemark® RP4.0 H1 CM LOC® FLEX abutment Brånemark® RP4.0 H2 CM LOC® FLEX abutment Brånemark® RP4.0 H3 CM LOC® FLEX abutment Brånemark® RP4.0 H4 CM LOC® FLEX abutment Brånemark® RP4.0 H5 CM LOC® FLEX case guide Brånemark® RP4.0 30°	05002114 05002115 05002116 05002117 05002118 05002188 05002261 05002262 05002263 05002264 05002265 05002309	35 Ncm -
Osstem® TS Regular	D 4.0, 4.5, 5.0, 6.0, 7.0	CM LOC® abutment TS Regular H1 CM LOC® abutment TS Regular H2 CM LOC® abutment TS Regular H3 CM LOC® abutment TS Regular H4 CM LOC® abutment TS Regular H5 CM LOC® case guide TS Regular 20°	05002159 05002160 05002161 05002162 05002163 05002197	30 Ncm -
Osstem® TS Mini	D 3.5	CM LOC® abutment TS Mini H1 CM LOC® abutment TS Mini H2 CM LOC® abutment TS Mini H3 CM LOC® abutment TS Mini H4 CM LOC® abutment TS Mini H5 CM LOC® case guide TS Mini 20°	05002590 05002591 05002592 05002593 05002594 05002596	30 Ncm
Osstem® US Regular	D 4.1	CM LOC® abutment US Regular H1 CM LOC® abutment US Regular H2 CM LOC® abutment US Regular H3 CM LOC® abutment US Regular H4 CM LOC® abutment US Regular H5 CM LOC® case guide US Regular 20°	05002154 05002155 05002156 05002157 05002158 05002196	30 Ncm -
Osstem® SS Regular	D 4.8	CM LOC® abutment SS Regular H1 CM LOC® abutment SS Regular H2 CM LOC® abutment SS Regular H3 CM LOC® abutment SS Regular H4 CM LOC® abutment SS Regular H5 CM LOC® case guide SS Regular 20°	05002149 05002150 05002151 05002152 05002153 05002195	30 Ncm -
Straumann® RN 4.8	RN D 4.8 mm	CM LOC® abutment Straumann® RN H1 CM LOC® abutment Straumann® RN H2 CM LOC® abutment Straumann® RN H3 CM LOC® abutment Straumann® RN H4 CM LOC® abutment Straumann® RN H5 CM LOC® case guide Straumann® RN 4.8 20° CM LOC® FLEX abutment Straumann® RN H1 CM LOC® FLEX abutment Straumann® RN H2 CM LOC® FLEX abutment Straumann® RN H3 CM LOC® FLEX abutment Straumann® RN H4 CM LOC® FLEX abutment Straumann® RN H5 CM LOC® FLEX case guide Straumann® RN 4.8 30°	05001408 05001409 05001410 05001411 05001412 05001578 05001720 02001721 05001722 05001723 05001724 05001884	35 Ncm -
Straumann® RC 4.1 / 4.8	RC D 4.1, 4.8 mm	CM LOC® abutment Straumann® RC H1 CM LOC® abutment Straumann® RC H2 CM LOC® abutment Straumann® RC H3 CM LOC® abutment Straumann® RC H4 CM LOC® abutment Straumann® RC H5 CM LOC® case guide Straumann® RC 20°	05001413 05001414 05001415 05001416 05001417 05001579	35 Ncm -

		CM LOC® FLEX abutment Straumann® RC H1 CM LOC® FLEX abutment Straumann® RC H2 CM LOC® FLEX abutment Straumann® RC H3 CM LOC® FLEX abutment Straumann® RC H4 CM LOC® FLEX abutment Straumann® RC H5 CM LOC® FLEX case guide Straumann® RC 30°	05001725 05001726 05001727 05001728 05001729 05001885	-
Straumann® NNC 3.5	NNC D 3.5 mm	CM LOC® abutment Straumann® NNC H1 CM LOC® abutment Straumann® NNC H2 CM LOC® abutment Straumann® NNC H3 CM LOC® abutment Straumann® NNC H4 CM LOC® abutment Straumann® NNC H5 CM LOC® case guide Straumann® NNC 20°	05001418 05001419 05001420 05001421 05001422 05001580	35 Ncm -
Sweden + Martina Premium Kohno	D 3.3 mm	CM LOC® abutment Premium Kohno 3.3 H1 CM LOC® abutment Premium Kohno 3.3 H2 CM LOC® abutment Premium Kohno 3.3 H3 CM LOC® abutment Premium Kohno 3.3 H4 CM LOC® abutment Premium Kohno 3.3 H5 CM LOC® case guide Premium Kohno 3.3 20°	05002119 05002120 05002121 05002122 05002123 05002189	30 Ncm -
	D 3.8 mm	CM LOC® abutment Premium Kohno 3.8 H1 CM LOC® abutment Premium Kohno 3.8 H2 CM LOC® abutment Premium Kohno 3.8 H3 CM LOC® abutment Premium Kohno 3.8 H4 CM LOC® abutment Premium Kohno 3.8 H5 CM LOC® case guide Premium Kohno 3.8 20°	05002124 05002125 05002126 05002127 05002128 05002190	-
	D 4.25 mm	CM LOC® FLEX abutment Premium Kohno 3.8 H1 CM LOC® FLEX abutment Premium Kohno 3.8 H2 CM LOC® FLEX abutment Premium Kohno 3.8 H3 CM LOC® FLEX abutment Premium Kohno 3.8 H4 CM LOC® FLEX abutment Premium Kohno 3.8 H5 CM LOC® FLEX case guide Premium Kohno 3.8 30°	05002266 05002267 05002268 05002269 05002270 05002310	-
	D 4.25 mm	CM LOC® abutment Premium Kohno 4.25 H1 CM LOC® abutment Premium Kohno 4.25 H2 CM LOC® abutment Premium Kohno 4.25 H3 CM LOC® abutment Premium Kohno 4.25 H4 CM LOC® abutment Premium Kohno 4.25 H5 CM LOC® case guide Premium Kohno 4.25, 5.0, 6.0 20°	05002129 05002130 05002131 05002132 05002133 05002192	-
	D 5.0, 6.0 mm	CM LOC® FLEX abutment Premium Kohno 4.25 H1 CM LOC® FLEX abutment Premium Kohno 4.25 H2 CM LOC® FLEX abutment Premium Kohno 4.25 H3 CM LOC® FLEX abutment Premium Kohno 4.25 H4 CM LOC® FLEX abutment Premium Kohno 4.25 H5	05002271 05002272 05002273 05002274 05002275	
	D 5.0, 6.0 mm	CM LOC® abutment Premium Kohno 5.0, 6.0 H1 CM LOC® abutment Premium Kohno 5.0, 6.0 H2 CM LOC® abutment Premium Kohno 5.0, 6.0 H3 CM LOC® abutment Premium Kohno 5.0, 6.0 H4 CM LOC® abutment Premium Kohno 5.0, 6.0 H5	05002134 05002135 05002136 05002137 05002138	
	D 5.0, 6.0 mm	CM LOC® FLEX abutment Premium Kohno 5.0, 6.0 H1 CM LOC® FLEX abutment Premium Kohno 5.0, 6.0 H2 CM LOC® FLEX abutment Premium Kohno 5.0, 6.0 H3 CM LOC® FLEX abutment Premium Kohno 5.0, 6.0 H4 CM LOC® FLEX abutment Premium Kohno 5.0, 6.0 H5 CM LOC® FLEX case guide Premium Kohno 4.25, 5.0, 6.0 30°	05002276 05002277 05002278 05002279 05002280 05002312	-
Zimmer MIS® (standard narrow) BioHorizons®	Tapered Screw-Vent® / MiS® Seven / BioHorizons® D 3.5 mm	CM LOC® abutment BioHorizons® 3.5 H1 CM LOC® abutment BioHorizons® 3.5 H2 CM LOC® abutment BioHorizons® 3.5 H3 CM LOC® abutment BioHorizons® 3.5 H4 CM LOC® abutment BioHorizons® 3.5 H5 CM LOC® case guide BioHorizons® 3.5 20°	05001462 05001463 05001464 05001465 05001466 05001589	30 Ncm -
	Tapered Screw-Vent® / MiS® Seven / BioHorizons® D 4.5 mm	CM LOC® abutment BioHorizons® 4.5 H1 CM LOC® abutment BioHorizons® 4.5 H2 CM LOC® abutment BioHorizons® 4.5 H3 CM LOC® abutment BioHorizons® 4.5 H4 CM LOC® abutment BioHorizons® 4.5 H5 CM LOC® case guide BioHorizons® 4.5 20°	05001467 05001468 05001469 05001470 05001471 05001590	-

		CM LOC® FLEX abutment BioHorizons® 4.5 H1	05001780	
		CM LOC® FLEX abutment BioHorizons® 4.5 H2	05001781	
		CM LOC® FLEX abutment BioHorizons® 4.5 H3	05001782	
		CM LOC® FLEX abutment BioHorizons® 4.5 H4	05001783	
		CM LOC® FLEX abutment BioHorizons® 4.5 H5	05001784	
		CM LOC® FLEX case guide BioHorizons® 4.5 30°	05001895	-

**Für Bestellungen und weitere Information steht Ihnen unser Kundenbüro gerne zur Verfügung.
For orders and further information, please contact the AUROSAN Customer Care Center.**

e-Mail service@aurosan.de | phone +49 201 219617-01 | fax +49 201 219617-31