

REVOIS®

Dental Implantat System

step ahead to future : REVOIS®

selected few parts - easy to use

ready for the digital age

best value in high quality

REVOIS® PRO



REVOIS® compact



At a glance: CAD/CAM, REVOIS® PRO & beyond

Digital Processing and CAD/CAM capabilities are today standard for a modern Dental Implant System.

For **REVOIS® PRO**, but also other implant systems, we and our partner **Cendres & Métaux** offer to you and your patients the the CAM-milling of highest-quality suprastructures at affordable pricing and with a 10-years guarantee.

For available 3D-data and information about compatible data formats, and to order scanbodies etc., please contact the AUROSAN Customer Care Center, or visit our Website at www.revois-dental.com.

The following types of suprastructures, compatible with **REVOIS®**, can be milled. All work is performed in Switzerland.

Indications	Abutments (lock)	Bridges, bars, full or partial arches; Wrap-arounds (no-lock)
Scanner / Software	Exocad	Immetrics
Materials	CrCo, Titanium; Gold upon request	
Time to delivery	1 - 2 working days	2- 3 working days

Efficient Working: You select – we deliver:

You scan and upload your STL-files onto www.cmsa.ch/cadcam-order, or you send us the impression model. In the latter, we scan and design, and send you the design-proposal for approval. Then we produce the work and mail it back to you.

Our guarantee lasts for 10 years for the milled suprastructure, and also the free-of-charge replacement in case of an acknowledged product defect. For details, please refer to www.cmsa.ch/cm-warranty.

Pricing (net): transparent and competitive

Indication	Material	Based on your scan data	Based on your model
Abutment	CoCo, Ti	120.00 €	160.00 €
bridge, directly screw-fixed	CoCo, Ti	145.00 €	185.00 €
Bar, directly screw-fixed	CoCo, Ti	145.00 €	185.00 €

Price is net and excludes transport. Ordering and invoicing via **AUROSAN**. Please contact us for any questions or further details.

AUROSAN GmbH

Customer Care Center
Frankenstr. 231
D-45134 Essen
Germany
Tel +49 201 219617-01
Fax +49 201 219617-31
eMail: service@aurosan.de



www.revois-dental.com



Cendres + Métaux SA

Scan+Design Center
Rue de Boujean 122
CH-2501 Biel-Bienne
Switzerland
Tel +41 58 360 23 20
Fax +41 58 360 20 10
eMail: dental-cadcam@cmsa.ch

www.cmsa.ch/cadcam-order

We look forward to your work order. Our quality and service will convince and satisfy you.

REVOIS®

Dental Implantat System

AUROSAN is sponsoring partner of the

ICOI EUROPE

INTERNATIONAL CONGRESS OF ORAL IMPLANTOLOGISTS